OSAMU-Retainer®

D

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The OSAMU-Retainer® has been developed by Dr. Osamu Yoshii, Tokyo, as treatment concept for the retention phase. After more than 2 years of clinical tests the retainer is now already in successful clinical application for many years. The retainer is fabricated individually for upper and lower jaw of a combination of highest quality pressure moulding materials, a soft BIOPLAST® foil (1.5 mm) and a hard-elastic IMPRELON® S pd foil (0.75 mm). The horse-shoe-shaped foils cover only the dental arch as well as parts of the mucosa in the apical area.

The tooth up to the alveolar area is covered by the soft BIOPLAST® foil while the occlusal surface is precisely held in the hard IMPRELON® S pd foil. These advantages in form and function of the OSAMU-Retainer® guarantee high retention results.

Other advantages are outstanding esthetics by the absolute transparency as well as the easy and efficient fabrication.

- The OSAMU-Retainer® guarantees maximum mechanical retention and eliminates the problem of relapse.
- The OSAMU-Retainer® is worn only at night. Nevertheless, speech with the splint is not impaired.
- It is easy for the patient to insert and remove the splint.
- Easy, low cost and time saving fabrication procedures.
- Slight movements and corrections of the teeth are possible – similar to the Positioner – if the OSAMU-Retainer® is fabricated on a Set-Up model.

Fabrication Procedure of an OSAMU-Retainer®

1. Draw approximate design of the Osamu-Retainer® on the plaster model and apply BIOPLAST® insulating agent to this area.
2. Heat the BIOPLAST® 1.5 mm foil and press it onto the model placed on the model platform.
3. Due to stability reasons cut off the BIOPLAST® foil only along the bottom of the cast by using the SD-foil scissors or the SD-scalpel.
4. Cut off the occlusal surfaces and incisal edges up to the contact points, finish the edges using the BIOPERM trimmer and clean the BIOPLAST® foil carefully with warm water.
5. Place cleaned BIOPLAST® foil back onto the model. For pressurizing IMPRELON® S pd (0.75 mm) embed the model in the pellets of the pressure moulding machine.
6. For an optimum compound of both foils apply Osamu-Bond on the complete BIOPLAST® foil during the last 20 sec. of the heating time of IMPRELON® S pd.
7. After the prolonged polymerization cut the OSAMU-Retainer® along the marked line.
8. Finishing and polishing of the marginal edges is done by using the Finishing Set/DIMO® Set.
9. Completed appliance, inserted in upper jaw.

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OSAMU-Retainer® - Hints for Fabrication Procedure

In order to pressure mould BIOPLAST® 1.5 mm we recommend to use only flat trimmed models to be placed on the model platform with the even side. For absolute transparency, we recommend the use of BIOPLAST® Insulating Agent.

For stability reasons, trim the pressed BIOPLAST® mould first of all to the bottom of the model. For trimming, you can either use SD-scissors A, SD-scalpel, diamond cutting disk or HM-cutting bur. Avoid any dirt on the BIOPLAST® sheet; before laminating BIOPLAST® with IMPRELON® S pd 0.75 mm make sure to clean the contact faces which

will be laminated together with lukewarm water and dry afterwards. For pressure moulding IMPRELON® S pd 0.75 mm embed the model in the pellets in order to reduce stretching and to guarantee the necessary stability of the sheet. Add OSAMU-Bond with a brush on the BIOPLAST® sheet during the last 10-15 seconds of heating time of IMPRELON® S pd 0.75 mm.

For long permanent compound, we recommend to prolong the cooling time of IMPRELON® S pd 0.75 mm and to cool the OSAMU-Retainer® 4-5 minutes under pressure.

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