

Ortho FlexTech®

Dental Retention System

INDICATIONS FOR USE: Intended for use as a Lingual Retainer Wire. This product is not indicated for individuals with sensitivity to Nickel.

INSTRUCTIONS FOR USE:

1. Prophy and sandblast the enamel, rinse and dry. Etch the enamel for 30 seconds. Rinse and Dry.
2. Apply one coat of Assure® Plus to the etched enamel. Lightly air dry.
3. Apply a minimal amount of L.C.R.™ GoTo™ or FlowTain™ to the lingual surface of each tooth.
 - a. Place the Ortho FlexTech® chain in the applied adhesive (“undercoat”) on the lingual surface. Position for best alignment and **take up any slack in the chain.**
 - b. Alternate Method: Place one end of Ortho FlexTech® chain in the applied adhesive paste (“undercoat”) on the lingual cuspid surface. Position for best alignment and light cure to hold chain in place. Proceed to the next tooth, **taking up any slack in the chain.** Repeat this process until the opposing cuspid has been reached.
4. Place an “overcoat” of L.C.R.™ GoTo™ or FlowTain™ paste on top of the cured material to cover the Ortho FlexTech™ chain. Light cure the “overcoat” paste on each bonded tooth.
5. Check for interference between upper and lower. Polish any rough spots or high points carefully.

NOTE: If used on upper to hold diastema, consider secondary .0175 twist 1 to 1 or Retainium™, as the chain can stretch slightly allowing space to reopen.



Do not Re-Use. This is a single use device

REORDER NUMBERS: REF# OFT2, Ortho FlexTech Gold 30”, REF# SSOFT, Ortho FlexTech Stainless Steel 30”, REF# OFT, Ortho FlexTech Gold Etched, 30”, REF# SSOFTE, Ortho FlexTech Stainless Steel Etched.

WARRANTY: Reliance Orthodontics Products, Inc. recognizes its responsibility to replace products if proven to be defective. Reliance Orthodontic Products, Inc. does not accept liability for any damages or loss, direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Rx Only: U.S. Federal law restricts this device to sale by or on the order of a dental professional.



RELIANCE[®]
ORTHODONTIC
PRODUCTS

1540 West Thorndale Ave., Itasca, IL 60143 U.S.A.
1-800-323-4348 · www.RelianceOrthodontics.com

MD

Scan for more information

EC

REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Australian Sponsor: Emergo
Australia
201 Sussex St. Darling Park
Tower II, Level 20,
Sydney, NSW 2000 Australia

Check out our new learning
platform Reliance U!
www.RelianceU.tech



WARNING: This product can expose you to Nickel (metallic), which is known to the State of California to cause cancer. For more information, go to www.P65Warnings.ca.gov