



IMPRELON® S

– THE solution for invisible splint retainer!

Our IMPRELON® S is the ideal solution to be used as a retention foil.


The improved material properties convince with extremely break- and abrasion-resistancy.

The foils with a thickness of 0.75 or 1.0 mm and a diameter of 125 mm are pre-dried and individually vacuum-packed for immediate use.

IMPRELON® S pd 125 mm


IMPRELON® S pd 0.75 x 125 mm REF 3322

IMPRELON® S pd 1.0 x 125 mm REF 3324




High and constant force

→ Keeps teeth in position



— IMPRELON® S
— DURAN®
— Competitor

IMPRELON® S provides a constant force with improved retention by holding the teeth longer in position. Here in comparison with DURAN® and a competitor's material.



— IMPRELON® S
— DURAN®
— Competitor

IMPRELON® S shows a high and constant restoring force during continuous stress and relief.

Improved material properties

→ Increased wear comfort

→ Low level of cracks and scratches

